

Automated VCMA & AHE Measurement System – AutoVA

AutoVA - Automated VCMA/(magnetoelectric coupling) - Anomalous Hall effect Measurement System.

The software is used to control :

1. The E-field applied for the VCMA/Magnetoelectric coupling. (with the ramping function).
 1. The E-fields can have multi-variants - (dwelling, pulser).
2. The electromagnetic field.
3. The Source Measure Unit (SMU) to apply one bias-current, and meanwhile measuring the Hall Voltage. Hysteresis loops can be obtained by looping the electromagnetic field.

The software has following features:

1. compatible with a large amount of SMUs, and instruments.
2. Experiments procedures visualised
3. GUI for graphing live data
4. managing queues of experiments.
5. Data auto-generation

Table of Contents:

1. GUI Introduction
 1. Block1

2. [Block2](#)

3. [Block3](#)

4. [Block4](#)


2. [Experiment Running & Data Example](#)

1. [Experiment Running Example](#)

2. [Daya Example](#)

3. [Controllable Instrument List](#)

GUI introduction



Block 1:

- **Maximum Voltage (mV):** This setting determines the strength of positive electric fields. Use a positive value when you want positive electric fields, and set it to 0 if you prefer negative fields.
- **Dwelling Time (s):** This is the duration for which the targeting electric fields will be applied and held.
- **Minimum Voltage (mV):** This setting controls the intensity of negative electric fields. Set a negative value when you want negative electric fields and 0 when you want positive fields.
- **Voltage Step (mV):** This is the increment used when ramping up the voltage to reach the target electric field strength.
- **Delay Time (ms):** This delay time is used during the voltage step when ramping up to the target electric field strength.
- **Compliance Current (uA):** This current is applied to protect the sample from high leakage currents.
- **Maximum Current Sent to Coil (A):** This current is sent to the coil to create looping electromagnetic (EM) fields for generating hysteresis loops.
- **Step (A):** This is the current increment used during the looping of EM fields.
- **Directory:** Specify the path where you want to save the data.

Block 2:

The sequencer is provided which allows users to queue a series of measurements with varying one, or more, of the parameters. This sequencer thereby provides a convenient way to scan through the parameter space of the measurement procedure.

The sequences can be extended and shortened using the buttons `Add root item`, `Add item`, and `Remove item`. The latter two either add an item as a child of the currently selected item or remove the selected item, respectively. To queue the entered sequence the button `Queue sequence` can be used. If an error occurs in evaluating the sequence text-boxes, this is mentioned in the logger, and nothing is queued.

```
- "Maximum Voltage", "[0]"
-- "Minimum Voltage", "[-1000]"
--- "Dwelling Time", "[540]"
-- "Minimum Voltage", "arange(-2000, -4500, -500)"
--- "Dwelling Time", "[60, 180, 300]"
```

Finally, it is possible to create a sequence file such that the user does not need to write the sequence again each time. The sequence file can be created by saving current sequence built within the GUI using the `Save sequence` button or directly writing a simple text file. Once created, the sequence can be loaded with the `Load sequence` button.

Block 3:

Graphical Display:

- live plotting for data
 - Dock1: IV curve when E-fields are applied
 - Dock 2: Coil currents (EM after calibration) looped as a function of the Hall Voltage measured
- Experiments log

Block 4:



A queue system for managing large numbers of experiments

Experiment Running & Data Example


Experiment Running Example

One experiment running example is shown as followings:


1. Parameters input:




2. Applying the E-field:




3. Looping the electromagnetic fields and receive the Hall voltages:



4. Applying another E-field and EM looping:



5. Making the sequencer (Automated Measurement):



Dataset Example

Two files will be generated: 1. the file containing all the data sets and one separate file containing the AHE data

	A	B	C	D	E	
1	#Procedure: < _main_ .IVProcedure>					
2	#Parameters:					
3	# Bias Current (A): 0.0004 A					
4	# Maximum current sent to coil (A): 4 A					
5	# Step (A): 0.2 A					
6	# Waiting Time (s): 0.1 s					
7	# Compliance Current: 1.05 uA					
8	# Delay Time: 100 ms					
9	# Devices: 1 No.					
10	# Maximum Voltage: 0 mV					
11	# Minimum Voltage: -1000 mV					
12	# Dwelling Time: 5 s					
13	# Voltage Step: 20 mV					
14	#Data:					
15	Voltage (V)	Current (A)	Resistance (ohm)	FieldCurrent	Hall Voltage (V)	
16	0	-3.99E-08	0	nan	nan	
17	-0.02	3.41E-08	-587237.85	nan	nan	
18	-0.04	-1.60E-08	2500048.44	nan	nan	
19	-0.06	1.67E-08	-3595973.2	nan	nan	
20	-0.08	2.68E-08	-2984820.7	nan	nan	
21	-0.1	-2.58E-08	3878960.9	nan	nan	
22	-0.12	2.75E-08	-4356306.7	nan	nan	
23	-0.14	8.95E-09	-15648557	nan	nan	
24	-0.16	-3.96E-08	4041897.3	nan	nan	
25	-0.18	3.85E-08	-4676591.6	nan	nan	
26	-0.2	-2.75E-08	7268522.1	nan	nan	
27	-0.22	2.81E-08	-7822241.7	nan	nan	
28	-0.24	-2.81E-08	8544507.27	nan	nan	
29	-0.26	3.31E-08	-7852610.1	nan	nan	
30	-0.28	-3.94E-08	7099494.08	nan	nan	
31	-0.3	2.79E-08	-10744516	nan	nan	
32	-0.32	-3.85E-08	8312059.66	nan	nan	
33	-0.34	4.01E-08	-8488506.6	nan	nan	
34	-0.36	-3.97E-08	9078302.63	nan	nan	
35	-0.38	2.90E-08	-13122589	nan	nan	
36	-0.4	-2.44E-08	16363739.8	nan	nan	
37	-0.42	2.08E-08	-20174442	nan	nan	

3	nan	nan	nan	0	0.5454546	
4	nan	nan	nan	0.19942	0.5575067	
5	nan	nan	nan	0.39949	0.5565896	
6	nan	nan	nan	0.59956	0.5569956	
7	nan	nan	nan	0.79963	0.5424007	
8	nan	nan	nan	0.9997	0.5587709	
9	nan	nan	nan	1.19912	0.5427202	
0	nan	nan	nan	1.3992	0.5608138	
1	nan	nan	nan	1.59927	0.558616	
2	nan	nan	nan	1.79934	0.5595382	
3	nan	nan	nan	1.99941	0.5612893	
4	nan	nan	nan	2.19948	0.5562355	
5	nan	nan	nan	2.39955	0.5590432	
6	nan	nan	nan	2.59962	0.5572236	
7	nan	nan	nan	2.79905	0.5545201	
8	nan	nan	nan	2.99912	0.5550638	
9	nan	nan	nan	3.19919	0.5549194	
0	nan	nan	nan	3.39926	0.5544061	
1	nan	nan	nan	3.59868	0.5558647	
2	nan	nan	nan	3.7994	0.5560151	
3	nan	nan	nan	3.99948	0.5578393	
4	nan	nan	nan	3.7994	0.5446938	
5	nan	nan	nan	3.59868	0.5395019	
6	nan	nan	nan	3.39926	0.5396025	
7	nan	nan	nan	3.19919	0.5450151	
8	nan	nan	nan	2.99912	0.5393209	
9	nan	nan	nan	2.79905	0.5540192	
0	nan	nan	nan	2.59898	0.5591146	
1	nan	nan	nan	2.39955	0.553203	
2	nan	nan	nan	2.19883	0.5532082	
3	nan	nan	nan	1.99941	0.5571428	

Controllable Instrument List

Controllable instrument is not limited as following, the equipments with VISA, GPIB adpater (most of the NI instruments) can be controlled.

Active Technologies

Active Technologies AWG-401x 1.2GS/s Arbitrary Waveform Generator

Advantest

Advantest R3767CG Vector Network Analyzer

Advantest R6245/R6246 DC Voltage/Current Sources/Monitors

